



Ref. : .....

الرقم : .....

Date : .....

التاريخ : .....

لجنة شراء لوازم والخدمات الاستشارية ( الفرعية )  
قرار الإحالة رقم (13/ 2022) للعطاء رقم (9/2022)  
المتضمن " توريد وتركيب وتشغيل أجهزة بحث علمي لكلية التمريض "

تاريخ : 2022/8/ 8

استناداً لنص المادة (35.أ) والمادة (36.أ.ب) من نظام المشتريات الحكومية رقم(8) لسنة 2022  
والمادة (29.ب.د.و) من التعليمات التنفيذية التابعة له وبناء على الرأي الفني قررت اللجنة ما يلي  
\*اعتماد الدراسة الفنية والمالية لغايات الإحالة.

12. إحالة البنود ذوات الأرقام (1,2,6,8,11,12,13,15,16) على السادة شركة الوراثة للتكنولوجيا الحيوية  
بمبلغ إجمالي وقدره (27,601.00) سبعة وعشرون ألفاً وستمئة وواحد دينار فقط لا غير، شامل الضريبة  
وحسب الكشف المرفق رقم (1).

13. إحالة البنود ذوات الأرقام (7 , 3+4) على مؤسسة برجوان للتجهيزات الطبية و المخبرية بمبلغ إجمالي  
وقدره (15056.80) خمسة عشر ألفاً وستة وخمسون ديناراً و 800—س فقط لا غير، شامل الضريبة  
وحسب الكشف المرفق رقم (2).

14. إحالة البند رقم (5) على شركة الجينوم الطبية بمبلغ إجمالي وقدره (5392.84) خمسة آلاف وثلاثمائة  
و اثنان وتسعون ديناراً و 840—ساً فقط لا غير، شامل الضريبة وحسب الكشف المرفق رقم (3)  
15. إحالة البنود ذوات الأرقام (9,10) على السادة المؤسسة الحديثة للتوريدات العلمية بمبلغ إجمالي  
وقدره (7,996) سبعة آلاف وتسعمائة وستة وتسعون ديناراً فقط لا غير، شاملاً الضريبة وحسب الكشف  
المرفق رقم (4).

\*\*\*بلغت قيمة قرار الإحالة الإجمالية مبلغ وقدره (56046.64) ستة وخمسون ألفاً وستة وأربعون ديناراً  
و 640—ساً فقط لا غير، شاملاً الضريبة.

- 1- تدفع الطوابع القانونية خلال عشرة أيام من تاريخ التبليغ بقرار الإحالة النهائي
- 2- تعتبر جميع المخاطبات السابقة واللاحقة والملاحق والتبليغات وقرارات اللجنة وعروض الشركة جزء لا يتجزأ من قرار الإحالة .
- 3 - تقدم كفالة حسن تنفيذ بقيمة (10%) من قيمة المواد المحالة .
- 4- يرجى الالتزام بالشروط المبينة في القرار.
- 5-تفرض غرامة تأخير قيمتها من (0.001 \_ 0.003) تصاعدياً حسب مدة التأخير استناداً لنص المادة (78/1) من تعليمات تنظيم إجراء المشتريات الحكومية لسنة (2022) وبحد أقصى (15%) من قيمة العقد/ الإحالة.
- 6-تقديم كفالة صيانة بنكية او شيك مصدق بنكي بقيمة (5%) من قيمة الأجهزة المحالة سارية لمدة سنة ميلادية من تاريخ التسليم.
- 7-يعتبر مبلغ الإحالة شاملاً لاي ضريبة او رسوم او جمارك .

كثف الاجهزة المحالة رقم (1) على شركة الوراثة للتكنولوجيا الحيوية بموجب قرار الاحالة رقم (13/2022)-الطعام رقم (9 / 2022)

رقم البند في الطعام	البند	مواصفات	الوحدة	الكمية المحالة	سعر الوحدة شامل الضريبة	السعر الاجمالي	المنشأ	البلد
1	<b>Cell counter Counting slide (2 boxes X 50 slides X 2 tests)</b>	<b>LUNAII™ Automated Cell Counter</b> Cat # L40002 Included with the instrument: 1. 2 box of counting slides (50 slides, 100 counts) 2. 1 pack of Trypan blue stain 3. 1 USB drive (including user manual) 4. Power cord and adapter (without printer) 5. Fully automated cell counting and cell viability analysis 6. Autofocus technology for fast and reliable cell counting 7. Various histograms and cell-size gating 8. Unmatched counting accuracy and speed 9. Built-in computer 10. The most affordable running cost currently on the market	Ea.	1	2445	2445	Logos Biosystems, Inc / Korea	الولايات متحدة
2	<b>Real time PCR</b>	Alt: 1 <b>PCR system</b> Cat#: Gentier 48E, 0 48 samples, 4 filters With PC pr Laptop Real-time PCR system Gentier 48E/48R Automatic PCR Analysis System is specially designed for small and medium-sized laboratories, mobile laboratories, field testing professionals and so on, to support absolute quantitative analysis, relative quantitative analysis, isothermal amplification, dissociation curve analysis, genotyping, and endpoint fluorescence analysis, etc. With its excellent performance and portability, it can make your experiments more easy-to-use, accurate and efficient.	Ea.	1	10000	10000	TIANLONG/ CHINA	الصين
6	<b>Single channel micropipette set</b>	A. Adjustable volume pipettes - Single channel AP10 -- 0.5-10 µl	Ea.	1	110	110	Biosan/Europe	الصين
		B. Adjustable volume pipettes - Single channel AP20 -- 2-20 µl	Ea.	1	110	110	Biosan/Europe	الصين
		C. Adjustable volume pipettes - Single channel AP50 -- 5-50 µl	Ea.	1	110	110	Biosan/Europe	الصين
		D. Adjustable volume pipettes - Single channel AP100 -- 10-100 µl	Ea.	1	110	110	Biosan/Europe	الصين
		E. Adjustable volume pipettes - Single channel AP200 -- 20-200 µl	Ea.	1	110	110	Biosan/Europe	الصين
		F. Adjustable volume pipettes - Single channel AP1000 -- 100-1000 µl	Ea.	1	110	110	Biosan/Europe	الصين
		G. 4-position stand	Ea.	1	30	30	Biosan/Europe	الصين



8	<p><b>High speed homogenizer</b></p>	<p>Omni general laboratory homogenizer (GLH 850) Cat# 20-020 General info Omni GLH General Laboratory Homogenizer The Omni GLH 850 General Laboratory Homogenizer is a powerful 850-watt homogenizer that was specifically designed to disrupt a wide range of samples in volumes ranging from 200 ul to 5 L. The Omni GLH is compatible with stainless steel generator probes ranging in diameter from 5 mm to 30 mm and 7 mm Omni Tips. The GLH quickly and efficiently processes a wide variety of sample and the 850-watt motor is capable of processing even the toughest samples. The GLH 850 features an advanced, integrated, push button digital control system that ensures accuracy and repeatability. The unit's closed loop variable speed controller (5,000 – 25,000 rpm) maintains selected speeds within +/- 100 rpm of the original set point. The digital timer can be set from seconds to minutes and includes a 30 second countdown timer, which is ideal for remote applications in fume hoods or biosafety cabinets. As soft start feature ensures smooth processing and minimizes splashing.</p> <p><b>A. 5 mm X 75 mm Flat Bottom (Fine) Generator Probe</b> Sku: G5-75 5 mm X 75 mm Flat Bottom (Fine) Generator Probe Stator design allows fine to medium samples to pass through the processing windows. Comes with PTFE upper and lower bearings. The entire probe can be autoclaved, or easily disassembled for cleaning. Optional Stainless steel upper bearings extend bearing life but are not autoclavable</p> <p><b>c. 10 mm X 115 mm Saw Tooth (Fine) Generator Probe</b> Sku: G10-115 10 mm X 115 mm Saw Tooth (Fine) Generator Probe Saw Teeth add extra cutting action and stator design allows fine to medium samples to pass through the processing windows. Comes with PTFE upper and lower bearings. The entire probe can be autoclaved, or easily disassembled for cleaning. Optional Stainless steel upper bearings extend bearing life but are not autoclavable</p>	Ea.	1	1083	1083	OMNI International /USA	حالا تقريباً
11	<p><b>Freeze dryer (lyophilizer)</b></p>	<p>Biobase Freeze Dryer (Table Top Type) BK-FDI 0PT  Features: 1. LCD touch screen. 2. With historical data inquiry function 3. With USB interface to export the data. 4. Large capacity ice condenser trap and without coil inside. 5. With cascade refrigeration technology. 6. Pre-freezing shelf can be used as guide barrel to speed up the drying rate. 7. Ice condenser trap and operation panel are made by stainless steel. 8. Drying chamber is transparent which is visual and safe. 9. Shelves are stainless steel, which can be adjustable per the requirements.</p>	Each	1	3550	3550	Biobase/China	حالا تقريباً

كشف الأجهزة المحالة رقم (1) على شركة الورثة للتكنولوجيا الحيوية بموجب قرار الاحالة رقم (13/2022)- للمطاء رقم (9 / 2022)

12	<b>Ice maker machine</b>	<p>FIM20 Ice maker                      Specifications:                      Item No. FIM20                      Ice making capacity(kg/24h) 20                      Ice storage capacity(kg) 10                      Ice shape Irregular small granular snowflake crushedice                      Input current(w) 180                      Ambient Temp. 10-25°C                      Water Temperature 10°C                      Power Supply AC220V±10%,50/60HZ; 110V±10%,60HZ                      Refrigeration type Air cooling                      Tank shell Stainless steel                      Compressor / refrigerant CFC-free/R134a                      Standard Accessories Water inlet pipe*1, Drainage pipe*1, Ice spoon*1, Sealing washer*2</p>	Ea.	1	800	800	Biobase/China	حالا عقوديا
----	--------------------------	--	-----	---	-----	-----	---------------	----------------

كثف الأجهزة المتاحة رقم (1) على شركة الوراثة للتكنولوجيا الحيوية بموجب قرار الإحالة رقم (13/2022)-الطباء رقم (9 / 2022)

13	<p><b>Water purification system</b></p> <p>Labagua Bio ultrapure water system cat# BS-070106-A06</p> <p>1. Labagua Bio - Point-of-use ultrafilter, UV sterilization module, TOC monitor Features:          1. Volumetric dispense - enables the user to set accurate dispensing volume for each dispense cycle. The dispense volume can be set either from the keyboard or by using "teaching" mode.          2. Water quality - embedded recirculation loop ensures stable premium water quality and enables practical elimination of Total Organic Carbon (TOC).          3. Low running costs - Performance of the deionization and polishing modules is constantly monitored. Monitoring algorithm enables cutting running costs, as replacement of the modules is requested only when service life is close to the end.          4. Total organic carbon (TOC) monitor - organic contaminants may not have effect on conductivity of water, so conductivity sensors cannot be used for TOC monitoring. Therefore, a special TOC monitoring module is needed to measure TOC level.          5. Color graphic LCD display - system component status is reflected on the display in an intuitive color pattern (Green/Yellow/Red).          6. System flowchart - shows all component status and water quality parameters at a glance.          7. The Labagua systems include:          1. Boost pump          2. Pre-filter set          3. Reverse osmosis module          4. Deionization module          5. Final stage polishing module          6. 30L storage tank with an integrated Grade 2 dispensing valve          7. Recirculation system</p>	Ea.	1	5500	5500	Biobase/China	حاليا تقريباً
15	<p><b>Digital caliber</b></p> <p>Digital Caliper Item #: 702660</p> <p>Description          Constructed of hardened stainless steel with a bright LCD display. Easy-access, push-button controls allows metric/English selection and resets readout to zero. Operates on a 1.5-V silver oxide battery (included) and designed with automatic shutoff.          Includes instructions and a sturdy case. Specifications: range, 0 to 150 mm (0 to 6"), resolution, 0.01 mm (0.0005"). Scale is marked in both millimeters and inch</p>	Ea.	1	75	75	Carolina Biologicals	حاليا تقريباً
16	<p><b>Microwave oven</b></p> <p>Microwave oven 35 - 45 L</p>	Each	1	200	200	--	حاليا تقريباً
			المجموع الإجمالي		27601		

كثف الاجهزة المحالة رقم (2) على مؤسسة برهان للتجهيزات الطبية والمخبرية لقرار الاحكام رقم (13/2022)-القطر رقم (9 / 2022)

no#	رقم البند في القطر	البند	مواصفات	الوحدة	الكمية المحالة	سعر الوحدة يشمل الضريبة	السعر الاجمالي	المنفذ	
1	3+4	3. Western blotting system 4. Gel Document Imaging System / Chemiluminescence, Western Blotting	<p><b>Cat No. WESTERFLOWLITE</b>                      The chemiLITE workflow package includes everything required to take a mini gel from sample loading to quantified ECL Image. The injection moulded omniPAGE Mini vertical gel tank runs up to 4 hand-cast gels and works perfectly with the wide range of pre-stained protein ladders.                      The powerPRO 3AMP power supply runs both the gel tank and the mini format semi-dry blotter.                      Also included in the package are reagents and accessories for blocking, including a blot roller, PVDF or nitrocellulose membrane, single-step blocking solution, and high sensitivity ECL substrate.                      All that is required are the antibodies for your target protein.                      Product includes:                      OmniPAGE Mini System with Caster,                      Blue Wide Range Protein Ladder (10-245kDa),                      Tris-Glycine SDS 10X – 1 liter,                      PowerPRO 3AMP Power Supply, 300V, 3000mA, 300W,                      Semi-Dry Mini, 10 x 10cm System,                      Roll PVDF 28 cm x 3 m, 0.22 µm,                      Blot-Absorbent Filter paper, 10 x 10cm, pack of 50,                      10cm acrylic roller for western blotting,                      EZEZ JUST BLOCK Western Blocking solution and signal enhancer, 500ml, Lumipro ECL substrate kit:                      50ml Luminol/enhancer solution;                      50ml Peroxide solution,                      chemiLITE Chemiluminescence Imaging System</p>	Ea	1	13,224.00	13,224.00	Cleaver/UK	
				PC from local market with 21" screen	Ea	1	696	696	Local
				OmniPETTE 8 Channel, 50-300µl	Ea	1	284.2	284.2	Cleaver/UK
				Multi-channel pipettes 8 channels 20-200µl	Ea	1	284.2	284.2	Cleaver/UK
2	7	Multichannel Micropipette set	Multi-channel pipettes 8 channels 0.5-10µl	Ea	1	284.2	284.2	Cleaver/UK	
			Multi-channel pipettes 8 channels 5-50µl	Ea	1	284.2	284.2	Cleaver/UK	
						<b>المجموع الاجمالي</b>	<b>15,056,80</b>		



كثف الأجهزة المحالة رقم (3) على شركة الجينوم الطبية لقرار الاحالة رقم (13/2022) - للعطاء رقم (9 / 2022)

رقم	المنفذ	السعر الاجمالي	سعر الوحدة	الكمية المحالة	الوحدة	مواصفات	البيد	رقم البيد في العطاء	no#
3-4	MicroDigital Co., Ltd. Korea	5392.84	5392.84	1	Each	Nabi – UV/VIS Nano Spectrophotometer can conveniently read samples in cuvette or microvolume, and also receive results from single wavelength and full spectrum Nabi is a standalone instrument which utilizes 7.0 inch touch-enabled L.C. All data can be transferred to PC, just with a USB. - Sample processing in cuvette and micro volume pedestal- Spectrometry analyzes single or full spectrum- portable size- maneuver/operation screen- back up with a USB Flash drive- Parameters: DNA/RNA, Protein, Cell, End point, Kinetic, Spectrum	<b>Nanodrop</b>	5	1
		المجموع الاجمالي	5392.84						

كثف الاجهزة المحالة رقم (4) على الموسسة الحيوية للتوريدات الطبية لقرار الاحالة رقم (13/2022) - صطاء رقم (9 / 2022)

رقم البند	رقم البند في النظام	وصف البند	مواصفات	الوحدة	الكمية المحالة	سعر الوحدة شامل الضريبة	المجموع الاجمالي	المغلف	مدة التسليم
1	9	Sonicator with probe (ultrasonic homogenizer)	500 Watt Ultrasonic Processor, 220V (includes 1/2" probe) VCXS500 Features: • Programmable • Variable power output • 10 hour timer • Pulse mode • Wattage and Energy display • Temperature monitoring • Energy sepoint Specifications Power Output 500 watt maximum Frequency 20KHz Ultrasonic Converter Diameter 2.5" (64mm)/Length 7.25" (183mm) Standard Probe Part Tip Diameter 0.5" (13mm) with replaceable tip Length 5.5" (139mm) Material Ti6Al4V titanium alloy Converter Cable Length 6 ft. (1.8m) Support Stand and converter clamp Alt 1: A Manual Lift Vertical condenser Rotary Evaporator EV400H-V	Each	1	4729	4729	Sonics /USA	3-6 months
2	10	Rotary evaporator	B. VP30 Vacuum Pump model Evaporation Flask 1000ml volume - 29/32 connection size included Receiving Flask 1000ml volume - 3/20 round connection size included C. Receiving Flask - Volume: 50ml - Joint size: 3/20 D. Receiving Flask - Volume: 250ml - Joint size: 3/20	Each	1	417	417	Labtech/ China	3-6 months
				Each	1	162	162	Labtech designed in Italy, assembled in China	3-6 months
					المجموع الاجمالي	7996			